

FreeLine™ Oxygen Mask

Medium Concentration, PVC-Free

Intended use:

Medium concentration masks are intended for administration of oxygen when precise control of inspired oxygen concentration is not mandatory. During inspiration the patient draws gases both from oxygen flowing into the mask through the tubing, as well as from room air via ports on the sides of the mask. The masks are connected directly to a compressed air or oxygen supply. The products are intended for single use only.



Product range:				
	REF*	Item	Size	Description
	93-103FL	FreeLine [™] Medium Concentration Oxygen Mask PVC-Free	Adult	Medium Concentration Oxygen Mask with soft cushion and curved nose seal, Headband, No tubing
	93-223FL	FreeLine [™] Medium Concentration Oxygen Mask PVC-Free	Pediatric	Medium Concentration Oxygen Mask with soft cushion and curved nose seal, Headband, No tubing

^{*}Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:
Contains latex? No
Contains phthalates? No
Contains bisphenols? No
Oxygen mask Polypropylene
Elastic strap Lycra
Soft sealingThermoplastic elastomer
Oxygen connector Polypropylene

Certifications:	
CE-marking? Yes (CE 1	639)
Product classification	(EU)
Quality system certified to	3485
Sterilization:	
Sterile?	. No
Chalf life 5 x	(OOYO

Packaging specifications:							
REF	Inner package		Outer package/Shipper				Pallet
	Qty, pcs	Dimensions, mm (LxW)	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt. kg	Volume m ³	Qty, pcs
93-103FL	1	180 x 220	30	550 x 350 x 170	1.34	0.033	720
93-223FL	1	140 x 200	30	400 x 350 x 140	0.94	0.020	1,260

EAN codes:				
REF	Inner package	Outer package		
93-103FL	(01)06945099321046	(01)46945099321044		
93-223FL	(01)06945099321039	(01)46945099321037		

Packaging material:			
Inner package	Poly bag - Low density polyethylene / Insert - 80g paper		
Outer package	Carton - Corrugated board / Sticker - Coated paper		



FreeLine[™] Oxygen Mask



Medium Concentration, PVC-Free

Storage and handling:		
<u> </u>	FRAGILE - Contents of the transport package are fragile and shall therefore be handled with care	
*	KEEP DRY	
<u> </u>	THIS WAY UP - Indicates correct upright position of the transport package	
8	DO NOT USE HAND HOOKS - Hooks are prohibited for handling the transport package	



Ningbo Shengyurui Medical Appliances Co., Ltd.

No. 138 Binhaisi Rd. Hangzhou Bay New Zone, 315336, Ningbo City Zhejiang Province, P.R. China

Made in China

EC REP

Shanghai Intl. Holding Corp. GmbH Eiffestrasse 80, 20537 Hamburg Distributed by



ConvaTec Limited First Avenue, Deeside Industrial Park Deeside, Flintshire, CH5 2NU, UK

www.convatec.com