

Intubation Stylets

Intended use:

Intubation Stylets are used as an aid to facilitate the insertion of the Tracheal Tube. The products are intended for single use only.



Product range:

REF*	CH size	O.D., mm (±0.2)	Length, mm (±10) (straight part)	Total length, mm (±10) (including bent part)
UM64200006	CH06	2.0	272	305
UM64200010	CH10	3.3	346	384
UM64200012	CH12	4.0	337	384
UM64200014	CH14	4.7	337	384

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
Contains phthalates? No
StyletAluminum
SheathPolyvinyl chloride

Certifications:

CE-marking? Yes (0123)
Product classification Class IIA (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN) 37469

Sterilization:

Sterile?Yes
Sterilized by Ethylene oxide
Shelf life 5 years

Packaging specifications:

REF	Primary package/ Single peelpack		Secondary package/ Shelf box		Tertiary package/ Shipper				Pallet Qty, pcs
	Qty, pcs	Dimensions, mm (LxW)	Qty, pcs	Dimensions, mm (LxWxD)	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt. kg	Volume m ³	
UM64200006	1	500 x 66	20	510 x 80 x 45	400	520 x 340 x 250	5.2	0.044	6,400
UM64200010	1	500 x 66	20	510 x 80 x 45	400	520 x 340 x 250	6.8	0.044	6,400
UM64200012	1	500 x 66	20	510 x 80 x 45	400	520 x 340 x 250	8.0	0.044	6,400
UM64200014	1	500 x 66	20	510 x 80 x 45	400	520 x 340 x 250	9.2	0.044	6,400





EAN codes:

REF	Primary package	Secondary package	Tertiary package
UM64200006	(01)30768455148470	(01)50768455148474	(01)70768455148478
UM64200010	(01)30768455148487	(01)50768455148481	(01)70768455148485
UM64200012	(01)30768455148685	(01)50768455148689	(01)70768455148683
UM64200014	(01)30768455148494	(01)50768455148498	(01)70768455148492

Packaging material:

Primary package	Peelback pouch - Coated paper, Polyethylene terephthalate/Polyethylene
Secondary package	Ivory cardboard
Tertiary package	Corrugated cardboard

Storage and handling:

	FRAGILE - Contents of the transport package are fragile and shall therefore be handled with care
	KEEP DRY
	THIS WAY UP - Indicates correct upright position of the transport package
	DO NOT USE HAND HOOKS - Hooks are prohibited for handling the transport package



Well Lead Medical Co., Ltd.

C-4 Jinhu Industrial Estate,
Hualong, Panyu, Guangzhou
511434 P.R.China

Made in China



Shanghai Intl. Holding Corp. GmbH

Eiffestrasse 80, 20537 Hamburg
Germany

Distributed by



ConvaTec Limited
First Avenue, Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

www.convatec.com

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled.
All trade marks and logos are the property of a ConvaTec group company.

© 2017 ConvaTec Inc. AP-017959-MM, ver. JUL 2017