


Frazier Suction Tube

Intended use:

Suction tubes are intended to be used with a suction handle for the removal of blood and fluids during surgery. The product is intended for single use only.



Product range:

	REF*	Description	Connector(s)	Wrapping	I.D., mm	O.D., mm	Tube length
	UM114.2.200	Frazier/Fine Instrument Sterile Suction Tubing	One white and one blue connector	Double wrapping	6	9	2m (6.5')
	UM114.2.300	Frazier/Fine Instrument Sterile Suction Tubing	One white and one blue connector	Double wrapping	6	9	3m (10')
	UM114.2.450	Frazier/Fine Instrument Sterile Suction Tubing	One white and one blue connector	Double wrapping	6	9	4.5m (15')

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
Contains phthalates? No
TubePolyvinyl chloride
ConnectorsPolyvinyl chloride

Certifications:

CE-marking? Yes (CE 0197)
Product classification Class I
Quality system certified toEN ISO 13485:2012

Sterilization:

Sterile? Yes
Sterilized by Ethylene oxide
Shelf life 5 years

Packaging specifications:

REF	Primary package/ Single peelpack		Secondary package/ Polybag		Tertiary package/ Shipper				Pallet
	Qty, pcs	Dimensions, mm (LxW)	Qty, pcs	Dimensions, mm (LxW)	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt. kg	Volume m ³	Qty, pcs
UM114.2.200	1	370 x 160	-	-	25	400 x 300 x 210	1.0	0.025	1,000
UM114.2.300	1	370 x 160	-	-	25	400 x 300 x 210	1.0	0.025	1,000
UM114.2.450	1	350 x 250	-	-	25	400 x 300 x 350	1.0	0.042	600





EAN codes:

REF	Primary package	Secondary package	Tertiary package
UM114.2.200	(01)30768455147688	-	(01)70768455147686
UM114.2.300	(01)30768455147695	-	(01)70768455147693
UM114.2.450	(01)30768455147701	-	(01)70768455147709

Packaging material:

Primary package	Sleeve - High density polyethylene, Peel pouch - Coated paper, Polypropylene/Polyethylene
Secondary package	Polyethylene bag - High density polyethylene
Tertiary package	Corrugated board

Storage and handling:

	THIS WAY UP - Indicates correct upright position of the transport package
	KEEP DRY
	FRAGILE - Contents of the transport package are fragile and shall therefore be handled with care
	AVOID DIRECT SUNLIGHT



GCMEDICA ENTERPRISE LTD.(WUXI)

Loujin Industrial Park, Shuofang
Wuxi, 214143, P.R. China

Made in China



Kingsmead Service Limited

19 Mezzanine Floor, 19-21 Crawford Street
London, W1H 1PJ, UK

Distributed by



ConvaTec Limited
First Avenue, Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

www.convatec.com

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled.
All trade marks and logos are the property of a ConvaTec group company.

© 2017 ConvaTec Inc. AP-017937-MM, ver. JUL 2017