

# Oxygen Mask


## Low Concentration, Vanilla Scented

**Intended use:**

Low concentration masks are intended for administration of oxygen when precise control of inspired oxygen concentration is not mandatory. During inspiration the patient draws gases both from oxygen flowing into the mask through the tubing, as well as from room air via ports on the sides of the mask. The masks are connected directly to a compressed air or oxygen supply. The products are intended for single use only.



**Product range:**

	REF	Item	Size	Description	Tube length
	93-101MMV	Low Concentration Oxygen Mask	Adult	Adult Low Concentration Oxygen Mask (elongated) with nose clip, No tubing	-

**Material specifications:**

**Contains latex?** . . . . . No  
**Contains phthalates?** . . . . . No  
**Contains bisphenols?** . . . . . No  
**Mask** . . . . . Polyvinyl chloride  
**Elastic strap** . . . . . Polyester  
**Nose clip** . . . . . Aluminum  
**Oxygen nipple connector** . . . . . High density polyethylene

**Certifications:**

**CE-marking?** . . . . . Yes (CE 0120)  
**Product classification** . . . . . IIA  
**Quality system certified to** . . . . . EN ISO 13485:2012

**Sterilization:**

**Sterile?** . . . . . No  
**Shelf life** . . . . . 5 years

**Packaging specifications:**

REF	Primary package	Secondary package/Polybag	Tertiary package/Shipper			Pallet
	Qty, pcs	Qty, pcs	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt., kg	Qty, pcs
93-101MMV	1	10	50	350 x 280 x 240	20.0	900

**EAN codes:**

REF	Primary package	Secondary package	Tertiary package
93-101MMV	(01)06945099313119	(01)26945099313113	(01)46945099313117

**Packaging material:**

Primary package	Poly Bag - Low density polyethylene / Insert - 80g Paper
Secondary package	Poly Bag - Low density polyethylene / Sticker - Coated Paper
Tertiary package	Sticker - Coated Paper / Carton - Corrugated Board

# Oxygen Mask

Low Concentration, Vanilla Scented



## Storage and handling:

	<b>FRAGILE</b> - Contents of the transport package are fragile and shall therefore be handled with care
	<b>KEEP AWAY FROM RAIN</b> - Transport package shall be kept away from rain
	<b>THIS WAY UP</b> - Indicates correct upright position of the transport package
	<b>DO NOT USE HAND HOOKS</b> - Hooks are prohibited for handling the transport package



**Ningbo Shengyurui Medical Appliances Co., Ltd.**

No. 138 Binhaisi Rd.  
Hangzhou Bay New Zone, 315336, Ningbo City  
Zhejiang Province, P.R. China

Made in China



**Shanghai Intl. Holding Corp. GmbH**

Eiffestrasse 80, 20537 Hamburg  
Germany

Distributed by



ConvaTec Limited  
First Avenue, Deeside Industrial Park  
Deeside, Flintshire, CH5 2NU, UK

[www.convatec.com](http://www.convatec.com)

**For further information on these products, please contact your local ConvaTec sales representative.**

With reservations for misprints and alterations in product specifications. Printed version not controlled.  
All trade marks and logos are the property of a ConvaTec group company.

© 2017 ConvaTec Inc. AP-017345-MM, ver. FEB 2017