


Duodenal Tubes

4 eyes

Intended use:

A gastro-duodenal tube is inserted through the nose or mouth via the oesophagus into the stomach or duodenum. It is used to drain gastric contents, decompress the stomach or duodenum, obtain a specimen of the gastric contents and administer medicine, liquids and/or enteral feeding. The product is intended for single use only.



List of product range:

Duodenal Tubes - 4 eyes		Levin (with x-ray line)			Ryles (with x-ray line & tip)
Tube Outer diameter, mm	CH size	X-ray line, open end	X-ray line, closed end	X-ray line, hard introducer/mandrin, closed end	X-ray tip/separate Luer adapter
		REF*	REF*	REF*	REF*
2.67	CH 08	1000118x	1001818x		1100118x
3.33	CH 10	1000218x	1001918x		1100218x
4.00	CH 12	1000318x	1002018x		1100318x
4.67	CH 14	1000418x	1002118x	1000818x	1100418x
5.33	CH 16	1000518x	1002218x	1000918x	1100518x
6.00	CH 18	1000618x	1002318x	1001018x	1100618x
6.67	CH 20	1000718x	1002418x		1100718x

All REF: Shore A hardness: 78, Effective length: 1210 mm

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
 Contains DEHP? No
 Tube Polyvinyl chloride
 Connector Polyvinyl chloride

Certification:

CE-marking? Yes (CE 0086)
 Product classification Class IIA (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN) Levin .38561
 Global Medical Device Nomenclature (GMDN) Ryles .38561

Sterilisation:

Sterile? Yes
 Sterilized by Ethylene oxide
 Shelf life 5 years

Duodenal Tubes

4 eyes

Packaging:

Tube Outer diameter mm	CH size	REF*	Shelf box		Shipping carton				MOQ pcs.
			Qty., pcs.	Dimensions, mm	Qty., pcs.	Dimensions, mm	Gross Wt. kg	Volume m ³	
LEVIN									
2.67	CH 08	1000118x	100	297 x 190 x 366	400	391 x 304 x 743	4.5	0.088	800
		1001818x	100	297 x 190 x 366	400	391 x 304 x 743	4.5	0.088	800
3.33	CH 10	1000218x	100	297 x 190 x 366	400	391 x 304 x 743	5.5	0.088	800
		1001918x	100	297 x 190 x 366	400	391 x 304 x 743	5.5	0.088	800
4.00	CH 12	1000318x	100	297 x 190 x 366	400	391 x 304 x 743	6.5	0.088	800
		1002018x	100	297 x 190 x 366	400	391 x 304 x 743	6.5	0.088	800
4.67	CH 14	1000418x	50	294 x 125 x 360	300	391 x 304 x 743	7.5	0.088	600
		1000818x	50	294 x 125 x 360	300	391 x 304 x 743	8.5	0.088	600
		1002118x	50	294 x 125 x 360	300	391 x 304 x 743	8.5	0.088	600
5.33	CH 16	1000518x	50	294 x 125 x 360	300	391 x 304 x 743	9.5	0.088	600
		1002218x	50	294 x 125 x 360	300	391 x 304 x 743	9.5	0.088	600
		1000918x	50	294 x 125 x 360	300	391 x 304 x 743	10.5	0.088	600
6.00	CH 18	1000618x	50	294 x 125 x 360	300	391 x 304 x 743	10.5	0.088	600
		1002318x	50	294 x 125 x 360	300	391 x 304 x 743	10.5	0.088	600
		1001018x	50	294 x 125 x 360	300	391 x 304 x 743	11.5	0.088	600
6.67	CH 20	1000718x	1	Single-packed	100	392 x 306 x 375	6.5	0.049	200
		1002418x	1	Single-packed	100	392 x 306 x 375	6.5	0.049	200
RYLES									
2.67	CH 08	1100118x	100	297 x 190 x 366	400	391 x 304 x 743	4.5	0.088	800
3.33	CH 10	1100218x	100	297 x 190 x 366	400	391 x 304 x 743	5.5	0.088	800
4.00	CH 12	1100318x	100	297 x 190 x 366	400	391 x 304 x 743	6.5	0.088	800
4.67	CH 14	1100418x	50	294 x 125 x 360	300	391 x 304 x 743	8.5	0.088	600
5.33	CH 16	1100518x	50	294 x 125 x 360	300	391 x 304 x 743	9.5	0.088	600
6.00	CH 18	1100618x	50	294 x 125 x 360	300	391 x 304 x 743	10.5	0.088	600
6.67	CH 20	1100718x	1	Single-packed	100	392 x 306 x 375	6.5	0.045	200

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.

