




Yankauer Suction Handle

Intended use:

A **Yankauer Suction Handle** connected to suction connecting tube is used to suction oropharyngeal secretions in order to prevent aspiration. A Yankauer can also be used to clear operative sites during surgical procedures and its suctioned volume counted as blood loss during surgery. The product is intended for single use only.



Product range:

	REF*	Description	Wrapping	Vent	Colour
	UM1802	Yankauer Suction Handle, Sterile	Double wrapping	-	Clear
	UM1801	Yankauer Suction Handle, Sterile	Single wrapping	-	Clear
	UM0800	Yankauer Suction Handle, Non-Sterile	Bulk packed	-	Clear
	UM1812	Yankauer Suction Handle, Sterile	Double wrapping	Yes	Blue
	UM1811	Yankauer Suction Handle, Sterile	Single wrapping	Yes	Blue

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
Contains phthalates? No
Tip K-Resin
Handle K-Resin

Certifications:

CE-marking? Yes (CE 0197)
Product classification Class IIA (EU)
Quality system certified to EN ISO 13485:2012

Sterilization:

Sterile? Yes
 No: UM0800
Sterilized by Ethylene oxide (Sterile codes only)
Shelf life 5 years

Packaging specifications:

REF	Primary package/ Single peelpack		Secondary package/ Polybag		Tertiary package/ Shipper				Pallet
	Qty, pcs	Dimensions, mm (LxW)	Qty, pcs	Dimensions, mm (LxW)	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt. kg	Volume m ³	Qty, pcs
UM1802	1	308 x 83	25	500 x 400	100	360 x 300 x 210	3.0	0.023	4,000
UM1801	1	308 x 83	25	500 x 400	100	360 x 300 x 210	3.0	0.023	4,000
UM0800	100*	490 x 350	-	-	300	400 x 300 x 210	5.0	0.025	12,000
UM1812	1	308 x 83	25	500 x 400	100	360 x 300 x 210	3.0	0.023	4,000
UM1811	1	308 x 83	25	500 x 400	100	360 x 300 x 210	3.0	0.023	4,000

* Not single packed, Bulk packed in a polybag.





EAN codes:

REF	Primary package	Secondary package	Tertiary package
UM1802	(01)30768455147565	(01)50768455147569	(01)70768455147563
UM1801	(01)30768455147572	(01)50768455147576	(01)70768455147570
UM0800	(01)30768455147589	-	(01)70768455147587
UM1812	(01)30768455147596	(01)50768455147590	(01)70768455147594
UM1811	(01)30768455147602	(01)50768455147606	(01)70768455147600

Packaging material:

Primary package	Sleeve - High density polyethylene (if double wrapped), Peel pouch - Coated paper, Polypropylene/Polyethylene
Secondary package	Polyethylene bag - High density polyethylene
Tertiary package	Corrugated board

Storage and handling:

	THIS WAY UP - Indicates correct upright position of the transport package
	KEEP DRY
	FRAGILE - Contents of the transport package are fragile and shall therefore be handled with care
	AVOID DIRECT SUNLIGHT



GCMEDICA ENTERPRISE LTD.(WUXI)

Loujin Industrial Park, Shuofang
Wuxi, 214143, P.R. China

Made in China

EC REP

Kingsmead Service Limited

19 Mezzanine Floor, 19-21 Crawford Street
London, W1H 1PJ, UK

Distributed by



ConvaTec

ConvaTec Limited
First Avenue, Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

www.convatec.com

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled.

All trade marks and logos are the property of a ConvaTec group company.

© 2017 ConvaTec Inc. AP-018350-MM, ver. OCT 2017