

100% Silicone Foley Catheter

Male, Tiemann, 2-Way, Single Packed


Intended use:

The 100% Silicone Foley Catheters are indicated for transurethral drainage of urine from the bladder. The adult 2-way Tiemann catheters can be retained in the human body for up to 90 days. The products are designed for male users and intended for single use only.



Product range:

Silicone Foley Catheter, Male, Tiemann, 2-Way, Single Packed

	REF*	Total Length mm (±10)	Balloon Capacity ml	CH/FG	Tip	Packaging
 <p>UM41161610S</p>	UM41161210S	420	10	12	Tiemann	1/30
	UM41161410S	420	10	14	Tiemann	1/30
	UM41161610S	420	10	16	Tiemann	1/30
	UM41161810S	420	10	18	Tiemann	1/30
	UM41162010S	420	10	20	Tiemann	1/30
	UM41162210S	420	10	22	Tiemann	1/30
	UM41162410S	420	10	24	Tiemann	1/30

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
Contains phthalates? No
Catheter shaft and funnel Silicone
X-ray line Silicone
Balloon Silicone
Tip Silicone, Silicone adhesive
Hard valve Polyvinyl chloride, Polypropylene, Acrylonitrile-butadiene-styrene, Silicone, Stainless steel
Sleeve Polypropylene

Certifications:

CE-marking? Yes (0123)
Product classification Class IIB (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN) 34917

Sterilization:

Sterile? Yes
Sterilized by Ethylene oxide
Shelf life 5 years

100% Silicone Foley Catheter

Male, Tiemann, 2-Way, Single Packed



Packaging specifications:

REF	Primary package/ Single peelpack		Secondary package/ Shelf box		Tertiary package/ Shipper				Pallet
	Qty, pcs	Dimensions, mm (LxW)	Qty, pcs	Dimensions, mm (LxWxD)	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt. kg	Volume m ³	Qty, pcs
UM41161210S	1	505 x 66	1	510 x 80 x 20	30	525 x 225 x 260	3.0	0.031	720
UM41161410S	1	505 x 66	1	510 x 80 x 20	30	525 x 225 x 260	3.0	0.031	720
UM41161610S	1	505 x 66	1	510 x 80 x 20	30	525 x 225 x 260	3.0	0.031	720
UM41161810S	1	505 x 66	1	510 x 80 x 20	30	525 x 225 x 260	3.0	0.031	720
UM41162010S	1	505 x 66	1	510 x 80 x 20	30	525 x 225 x 260	3.0	0.031	720
UM41162210S	1	505 x 66	1	510 x 80 x 20	30	525 x 225 x 260	3.0	0.031	720
UM41162410S	1	505 x 66	1	510 x 80 x 20	30	525 x 225 x 260	3.0	0.031	720





EAN codes:

REF	Primary package	Secondary package	Tertiary package
UM41161210S	(01)30768455153931	(01)50768455153935	(01)70768455153939
UM41161410S	(01)30768455153948	(01)50768455153942	(01)70768455153946
UM41161610S	(01)30768455153955	(01)50768455153959	(01)70768455153953
UM41161810S	(01)30768455153962	(01)50768455153966	(01)70768455153960
UM41162010S	(01)30768455153979	(01)50768455153973	(01)70768455153977
UM41162210S	(01)30768455153986	(01)50768455153980	(01)70768455153984
UM41162410S	(01)30768455153993	(01)50768455153997	(01)70768455153991

Packaging material:

Primary package	Sleeve - Ethylene-vinyl acetate, Peelpack pouch - Coated paper, Polyethylene terephthalate/Polyethylene
Secondary package	Ivory cardboard
Tertiary package	Corrugated board

Storage and handling:

	FRAGILE - Contents of the transport package are fragile and shall therefore be handled with care
	KEEP DRY
	THIS WAY UP - Indicates correct upright position of the transport package
	DO NOT USE HAND HOOKS - Hooks are prohibited for handling the transport package


Well Lead Medical Co., Ltd.
 C-4 Jinhu Industrial Estate
 Hualong, Panyu, Guangzhou
 511434 P.R. China

Made in China

EC REP

Shanghai Intl. Holding Corp. GmbH
 Eiffestrasse 80, 20537 Hamburg
 Germany

Distributed by



ConvaTec Limited
 First Avenue, Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

www.convatec.com

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled.
 All trade marks and logos are the property of a ConvaTec group company.

© 2018 ConvaTec Inc. AP-019442-MM, ver. OCT 2018