

# Careline™ 100% Silicone Foley Catheter

## Male, Tiemann, 2-Way, with Prefilled Syringe


### Intended use:

The 100% Silicone Foley Catheters are indicated for transurethral drainage of urine from the bladder. The adult 2-way Tiemann catheters can be retained in the human body for up to 90 days. The products are designed for male users and intended for single use only.




### Product range:

#### Careline™ Silicone Foley Catheter, Male, Tiemann, 2-Way, Prefilled Syringe (Sterile 10% Aqueous Glycerine Solution), Includes Additional Emptying Syringe

	REF*	Total Length, mm (±10)	Balloon Capacity, ml	CH/FG	Tip	Packaging
 UM41161610P	UM41161210P	420	10	12	Tiemann	5/250
	UM41161410P	420	10	14	Tiemann	5/250
	UM41161610P	420	10	16	Tiemann	5/250
	UM41161810P	420	10	18	Tiemann	5/250
	UM41162010P	420	10	20	Tiemann	5/250
	UM41162210P	420	10	22	Tiemann	5/250
	UM41162410P	420	10	24	Tiemann	5/250

#### Careline™ Silicone Foley Catheter, Male, Tiemann, 2-Way, Prefilled Syringe (Sterile 10% Aqueous Glycerine Solution), Includes Additional Emptying Syringe, Single Packed

 UM41161610PS	UM41161210PS	420	10	12	Tiemann	1/30
	UM41161410PS	420	10	14	Tiemann	1/30
	UM41161610PS	420	10	16	Tiemann	1/30
	UM41161810PS	420	10	18	Tiemann	1/30
	UM41162010PS	420	10	20	Tiemann	1/30
	UM41162210PS	420	10	22	Tiemann	1/30
	UM41162410PS	420	10	24	Tiemann	1/30

\* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Distributed by


**ConvaTec**

For further information on these products, please contact your local ConvaTec sales representative.

# Careline™ 100% Silicone Foley Catheter

## Male, Tiemann, 2-Way, with Prefilled Syringe



### Material specifications:

**Contains latex?** . . . . . No  
**Contains phthalates?** . . . . . No  
**Catheter shaft and funnel.** . . . . . Silicone  
**Balloon** . . . . . Silicone  
**X-ray line** . . . . . Silicone  
**Tip** . . . . . Silicone, Silicone adhesive  
**Hard valve** . . . . . Polyvinyl chloride, Polypropylene,  
 Acrylonitrile-butadiene-styrene, Silicone, Stainless steel  
**Sleeve** . . . . . Polypropylene  
**Syringe** . . . . . Polypropylene  
**Glycerine solution** . . . . . Glycerin, Water

### Certifications:

**CE-marking?** . . . . . Yes (0123)  
**Product classification** . . . . . Class IIB (EU)  
**Quality system certified to** . . . . . ISO 13485  
**Global Medical Device Nomenclature (GMDN)** . . . . . 34917

### Sterilization:

**Sterile?** . . . . . Yes  
**Sterilized by** . . . . . Ethylene oxide  
**Shelf life** . . . . . 3 years

### Packaging specifications:

REF	Primary package/ Single peelpack		Secondary package/ Shelf box		Tertiary package/ Shipper				Pallet
	Qty, pcs	Dimensions, mm (LxW)	Qty, pcs	Dimensions, mm (LxWxD)	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt. kg	Volume m <sup>3</sup>	Qty, pcs
UM41161210P	1	505 x 80	5	510 x 80 x 75	250	775 x 525 x 430	9.0	0.175	1,000
UM41161410P	1	505 x 80	5	510 x 80 x 75	250	775 x 525 x 430	9.0	0.175	1,000
UM41161610P	1	505 x 80	5	510 x 80 x 75	250	775 x 525 x 430	9.0	0.175	1,000
UM41161810P	1	505 x 80	5	510 x 80 x 75	250	775 x 525 x 430	9.0	0.175	1,000
UM41162010P	1	505 x 80	5	510 x 80 x 75	250	775 x 525 x 430	9.0	0.175	1,000
UM41162210P	1	505 x 80	5	510 x 80 x 75	250	775 x 525 x 430	9.0	0.175	1,000
UM41162410P	1	505 x 80	5	510 x 80 x 75	250	775 x 525 x 430	9.0	0.175	1,000
UM41161210PS	1	505 x 80	1	510 x 80 x 25	30	525 x 275 x 260	3.0	0.038	720
UM41161410PS	1	505 x 80	1	510 x 80 x 25	30	525 x 275 x 260	3.0	0.038	720
UM41161610PS	1	505 x 80	1	510 x 80 x 25	30	525 x 275 x 260	3.0	0.038	720
UM41161810PS	1	505 x 80	1	510 x 80 x 25	30	525 x 275 x 260	3.0	0.038	720
UM41162010PS	1	505 x 80	1	510 x 80 x 25	30	525 x 275 x 260	3.0	0.038	720
UM41162210PS	1	505 x 80	1	510 x 80 x 25	30	525 x 275 x 260	3.0	0.038	720
UM41162410PS	1	505 x 80	1	510 x 80 x 25	30	525 x 275 x 260	3.0	0.038	720

### EAN codes:

REF	Primary package	Secondary package	Tertiary package
UM41161210P	(01)30768455150404	(01)50768455150408	(01)70768455150402
UM41161410P	(01)30768455150411	(01)50768455150415	(01)70768455150419
UM41161610P	(01)30768455150428	(01)50768455150422	(01)70768455150426
UM41161810P	(01)30768455150435	(01)50768455150439	(01)70768455150433
UM41162010P	(01)30768455150442	(01)50768455150446	(01)70768455150440
UM41162210P	(01)30768455150459	(01)50768455150453	(01)70768455150457
UM41162410P	(01)30768455150466	(01)50768455150460	(01)70768455150464
UM41161210PS	(01)30768455150473	(01)50768455150477	(01)70768455150471
UM41161410PS	(01)30768455150480	(01)50768455150484	(01)70768455150488
UM41161610PS	(01)30768455150497	(01)50768455150491	(01)70768455150495
UM41161810PS	(01)30768455150503	(01)50768455150507	(01)70768455150501
UM41162010PS	(01)30768455150510	(01)50768455150514	(01)70768455150518
UM41162210PS	(01)30768455150527	(01)50768455150521	(01)70768455150525
UM41162410PS	(01)30768455150534	(01)50768455150538	(01)70768455150532

Distributed by



For further information on these products, please contact your local ConvaTec sales representative.





# Careline™ 100% Silicone Foley Catheter Male, Tiemann, 2-Way, with Prefilled Syringe



## Packaging material:

Primary package	Sleeve - Ethylene-vinyl acetate, Peelpack pouch - Coated paper, Polyethylene terephthalate/Polyethylene
Secondary package	Ivory cardboard
Tertiary package	Corrugated board

## Storage and handling:

	<b>FRAGILE</b> - Contents of the transport package are fragile and shall therefore be handled with care
	<b>KEEP DRY</b>
	<b>THIS WAY UP</b> - Indicates correct upright position of the transport package
	<b>DO NOT USE HAND HOOKS</b> - Hooks are prohibited for handling the transport package



**Well Lead Medical Co., Ltd.**

C-4 Jinhu Industrial Estate  
Hualong, Panyu, Guangzhou  
511434 P.R. China

Made in China



**Shanghai International Holding  
Corp. GmbH (Europe)**

Eiffestrasse 80, 20537 Hamburg  
Germany

Distributed by



ConvaTec Limited  
First Avenue, Deeside Industrial Park  
Deeside, Flintshire, CH5 2NU, UK  
[critcare.convatec.com](http://critcare.convatec.com)

**For further information on these products, please contact your local ConvaTec sales representative.**

With reservations for misprints and alterations in product specifications. Printed version not controlled.

© 2021 ConvaTec. ®/™ indicate trademarks of the ConvaTec group of companies.

© 2021 ConvaTec Inc. AP-033279-EM-NP-SP, ver. DEC 2021